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1. **Description**
The Optimos™ Guidewire is used for the insertion of endoscopic devices, especially for use in selective cannulation of the biliary ducts including, but not limited to, the common bile, cystic, right and left hepatic ducts and for interventional EUS (endoscopic ultrasound) procedures of the gastrointestinal tract. The guidewire consists of a nitinol core with a PTFE jacket to provide endoscopic visualization of wire movement. A radiopaque TPU (thermoplastic polyurethane) loaded jacket encapsulates the distal tip of the guidewire core for radiopacity. The thermoplastic polyurethane loaded jacket of the guidewire is then coated with a hydrophilic coating.
2. **Indications for Use**
Used to assist in the insertion of endoscopic devices for endoscopy and ultrasonography
- **Patient Population**
Not specialized for a particular population, but applicable to all patients over 18 years of age; however, the decision rests with the physicians
3. **Performance and Clinical Benefits**
 - Performance : Biliary cannulation
 - Clinical benefit : Drainage success
4. **Contraindications**
This product is not to be used for interventional procedures, so do not use it for vascular procedures.
5. **Potential Complications**
 - Perforation
 - Bleeding
 - Pancreatitis
 - Bile peritonitis

* Potential complication related to ERCP, such as allergic reactions to the contrast agent.
6. **Precautions**
 - This product is designed for single use only so it must be neither re-sterilized nor reused.
 - This product is only to be used when it has not expired.
 - Procedures using this product must be performed only by appropriately trained doctors.
 - Check that other devices for the procedure are suitable to the guidewire (use of metal and/or sharp tips may result in damage to the external coating tip of the guidewire).
 - Ensure the diameter of the guidewire on the label and select the endoscopic devices with proper inner sizes.
 - Move the guidewire gently when passing it through the accessory channel of the endoscope.
 - Do not rub the tip of the guidewire with dry fabric.

경고 If any serious incident occurs in relation to this device, the user and/or patient should report to the competent authority of the Member State in which the user and/or patient is established and to the manufacturer.

7. **Warnings**
 - If the wire gets caught in strictures, push or twist the wire gently.
 - Keeping the wire in an unusual shape has a bad effect on its appearance.
 - The hydrophilic coating can be removed if an alcohol-based agent is used for the procedure.
 - This product is disposable and should therefore not be reused. If so, risks such as infection may occur.
 - The hydrophilic coating can be removed if this product is reused. This deterioration of hydrophilic function could pose a risk to a patient.
 - Do not re-sterilize or reuse this product.
 - This product is supplied in sterilized packages; do not use this product if its internal package has been damaged due to impact, etc.
 - Do not use this product when the patient has a severe fever, hemorrhaging, or difficulty in breathing.

8. **Instructions for use**
 - (1) Before you remove the guidewire from the dispenser, connect the syringe to the dispenser and flush with sterile water to soak the wire.
 - (2) Flush the accessory channel of the endoscope and the devices used for the procedure.

[ERCP Procedure]

 - (3) Insert the guidewire into the duct while monitoring fluoroscopically.
 - (4) Advance the device through the guidewire after confirming that the wire is located at the right spot.

[EUS Procedure]

 - (3) Advance the EUS FNA needle to puncture the target after confirming that the endoscope is located at the right spot and inject the contrast agent.
 - (4) Advance the guidewire through the FNA needle and pass the endoscopic device through the guidewire to the target lesion.
 - (5) Dispose of the wire after the procedure in accordance with the prescribed disposal procedure.
9. **Handling / Storage**
The expiration period for this product is three years. Do not use any part of this product if it has expired. Store this product in a dry, clean place at the temperature of 15°C~30°C. Keep the product away from direct sunlight.

Reuse Precaution Statement
This product has been sterilized using ethylene oxide (E.O.) gas. Do not use if the sterile package is opened or damaged. The company is not responsible for re-sterilized products.

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