UVENTA™

Urethral Stent

Table of Contents

English.....	3
Français.....	7
Deutsch	11
Italiano.....	15
Español	19
Português	23
Türkçe.....	27
Čeština.....	31
Dansk.....	35
Ελληνικά.....	39
Lietuvių.....	43
Polski.....	47
Română.....	51
Svenska.....	55
Norsk.....	59

MRI Information

MR Conditional

UVENTA™ Urethral Stent was established to be MR-conditional.

Non-clinical testing demonstrated that Uventa™ Urethral Stent is MR Conditional according to ASTM F2503. A patient with this device can be scanned safely, immediately after placement under the following Conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less

- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, UVENTA™ Urethral Stent produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, and horizontal field scanner) and 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

	<u>1.5-Tesla</u>	<u>3-Tesla</u>
MR system reported, whole body averaged SAR	2.9-W/kg	2.9-W/kg
Calorimetry measured values, whole body averaged SAR	2.1-W/kg	2.7-W/kg
Highest temperature change	+1.9°C	+2.4°C

These temperature changes will not pose a hazard to a human subject under the conditions indicated above.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of Uventa™ Urethral Stent. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 10-mm relative to the size and shape of this implant.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	7,695-mm ²	676-mm ²	7,984-mm ²	977-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

User's Manual

1. Description

Uventa™ Urethral Stent consists of the implantable full covered metallic Stent and introducer system. Uventa™ Urethral Stent is made of Nitinol wire. It is a flexible, fine mesh tubular prosthesis which has radiopaque markers on each end and at the center.

Model Name
Uventa™ Urethral Stent

The Stent is loaded in the introducer system. Upon deployment, the Stent imparts an outward radial force on the luminal surface of the spongy urethra to establish patency (See Figure 1).


Figure 1. Introducer System

2. Operation

The outer sheath is pulled back by immobilizing the hub in one hand, grasping the Y-connector with the other hand, and gently sliding the Y-connector along the pusher towards the hub. Retraction of the outer sheath releases the Stent.

3. Purpose of Use

Uventa™ Urethral Stent is intended for maintaining urethral patency in benign and malignant urethral obstruction.

WARRANTY

Taewoong Medical Co., Ltd. warrants that reasonable care has been applied within the design and subsequent manufacturing process of this instrument was applied. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by the operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storing, cleaning and sterilizing of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Taewoong's control directly affect the instrument and the results obtained from its use. Taewoong's obligation under this warranty is limited to the replacement of this instrument and Taewoong shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. Taewoong neither assumes, nor authorizes any other person to assume for it, and any other or additional liability or responsibility in connection with this instrument. Taewoong assumes no liability with respect to instruments reused, reprocessed or re-sterilized, and makes no warranties, expressed or implied, including but not limited to, merchantability or fitness for a particular purpose, with respect to such instruments.

4. Contraindication

The use of Uventa™ Urethral Stent is prohibited for, including but not limited to:

- Patients who had radical retropubic prostatectomy due to prostate cancer
- Patients with serious neurogenic bladder
- Patients with septicemia due to urinary tract infection
- Patients who showed adverse metallic Stent reaction (Serious urinary tract infection, repeatable migration of Stent, or foreign body sensation) after installing Stent
- Patients with urethral perforation or abscess (Exceptionally, metallic Stent can be installed for treatment of urethral perforation by the doctor's discretion)
- Patients with suspected or impending perforation
- Patients who have been treated due to acute urethral infection
- Recapturing a stent during its deployment is contraindicated.

5. Warnings

- The Stent contains nickel, which may cause allergic reactions in individuals with nickel sensitivity.
- Do not expose the introducer system to organic solvent (e.g. Alcohol)
- Uventa™ Urethral Stent can be removed when the doctor's clinical assessment indicates that the Stent can be safely removed. Caution should be exercised in deciding to and when removing the Stent.
- Uventa™ Urethral Stent cannot be removed when tumor in-growth/over-growth/occlusion of the Stent lumen is present.
- Uventa™ Urethral Stent can be repositioned immediately after deployment.
- Caution should be exercised not to move the position of the Stent when checking the position of the Stent with cystoscope after surgery
- Sexual activity should be prohibited while Stent is installed.
- When selecting the position of the Stent, be careful not to install the Stent in the external sphincter.
- Patients who use erection assistant devices or are prescribed injections for erectile

dysfunction should not use the Stent.

- Patients who were implanted with artificial urinary sphincter should not use the Stent.
- Do not attempt to recapture/reload a stent once its deployment is advanced.

6. Precautions

Read the entire User's Manual thoroughly before using this device. It should only be used by or under the supervision of physicians thoroughly trained in the placement of Stents. A thorough understanding of the techniques, principles, clinical applications and risks associated with this procedure is necessary before using the device.

- Care should be taken when removing the introducer system and guidewire immediately after Stent deployment since this may result in Stent dislodgement if the Stent has not been adequately deployed.
- Care should be taken when performing dilation after the Stent has been deployed as this may result in perforation, bleeding, Stent dislodgement or Stent migration.
- Expiration date indicated on the packing should be checked and the use of the product which had gone beyond its expiration date is prohibited.
- If the packaging is open or damaged, you cannot use it as Uventa™ Urethral Stent is supplied sterilized.
- This product is disposable medical equipment. Reuse of it after resterilization is prohibited.

7. Instructions for Damages

WARNING: Visually inspect the system for any sign of damage. DO NOT USE the system if any visible signs of damage are present. Failure to observe this precaution may result in patient injury.

8. Potential complications

Potential complications associated with the use of Uventa™ Urethral Stent may include, but is not limited to:

Procedural Complications

- Bleeding
- Pain
- Stent misplacement
- Inadequate expansion of Stent
- Urethral perforation

Post Stent Placement Complications

- Bleeding
- Pain
- Urethral fistula or perforation
- Incontinence
- Stent misplacement or migration
- Stent occlusion due to encrustation
- Stent occlusion due to tissue ingrowth
- Stent occlusion due to tissue overgrowth or hyperplasia
- Fever
- Foreign body sensation
- Sepsis
- Infection
- Tissue erosion around urethra due to Stent stimulation

9. Required equipments

- Urethroscope or Cystoscope (over 10Fr. Working Channel)
- Sizing catheter
- Fluoroscope / C-Arm (if necessary)
- Forceps (if necessary)
- 0.035" / 0.89 mm Guide wire

10. Procedure

① Examine Stricture

- a) Carefully examine both proximal and distal segment of stricture by endoscopy and/or fluoroscopy.
- b) The internal luminal diameter should be measured accurately with endoscope and/or fluoroscopy.

② Determine Stent Size

- a) Measure the length of target stricture.
- b) Select a Stent size that is longer than the measured length of target stricture in order to fully cover both ends of the lesion (more than 5mm in the direction of the sphincter, more than 10mm in the direction of the external urethral orifice).

③ Stent Deployment Preparation

- Uventa™ Urethral Stent can be placed with the aid of fluoroscopy and/or endoscopy.
 - Pass a 0.035" (0.89 mm) guide wire to the region of the stricture.
- a) Pass the guide wire through the place where Stent will be positioned with the aid of fluoroscopy and/or endoscopy.
 - b) Ensure that the valve of Y-connector connecting the inner sheath and outer sheath is locked by rotating proximal valve end in a clockwise direction to prevent premature Stent deployment.

- c) Flush the inner lumen of introducer system with the sterilized saline solution.

④ Stent Deployment Procedure


Figure 2

- a) Under the fluoroscope and/or endoscopic guidance, position the introducer system to the center of the target stricture.
- b) Once the introducer system is in the correct position for deployment, unlock the proximal valve of the Y-connector by turning the valve more than twice in the anti-clockwise direction.
- c) To begin Stent deployment, immobilize the hub in one hand and grasp the Y-connector with the other. Gently slide the Y-connector back along the pusher towards the hub.
- d) When the center X-ray marker reaches the center of the target stricture, continue pulling the Y-connector until the Stent is fully deployed. (See figure 3)


Figure 3

CAUTION Do not push or pull on the hub with the Stent partially deployed. The hub must be securely immobilized. Inadvertent movement of the hub may cause misalignment of the Stent and possible damage to urethral tissue.

⑤ After Stent Deployment

- a) To confirm expansion, examine the Stent by fluoroscopy and/or endoscopy.
- b) Carefully remove the introducer system, guide wire and endoscope from the patient. If excessive resistance is felt during removal, wait for 3~5 minutes to allow further Stent expansion. (Place the inner sheath back into the outer sheath as the original state prior to removal.)
- c) Balloon dilatation inside the Stent can be performed with the discretion of the doctor.
- d) Wash with about 50cc of warm saline solution (38~45°C) after surgery to maximize Stent expansion.

11. Perform post-implant routine procedures

- a) Assess the size and stricture of the Stent lumen. A Stent may require 1 to 3 days to expand fully.
- b) Doctor's experience and discretion can determine the appropriate drug or procedure for each patient.
- c) After implantation, patient should follow the prescription of the treating doctor.
- d) Observe the patient for development of any complications.

12. Instructions for relocation and removal of Uventa™ Urethral Stent

Visually examine the Stent for any tumor in-growth/over-growth into the Stent lumen or whether the Stent is occluded. If the Stent lumen is clear, carefully remove the Stent by using forceps and/or an appropriate device. When removing, cleansing the Stent with cold saline solution (less than 4°C) can make the removal easier.

Grasp the retrieval string and/or collapse the proximal end of the Stent, and then carefully retrieve the Stent. If the Stent cannot be easily withdrawn, stop removing the Stent and make inquiry to the manufacturer or seller.

Caution: Do not use excessive force to remove the Stent as it may cause disconnect to the retrieval string.

To relocate Uventa™ Urethral Stent immediately after deployment, use forceps or an appropriate device to grasp the retrieval string and gently adjust to the correct placement. Adjustment should be made only toward the external urethral orifice.

The timing of the removal should be determined by the doctor at his/her discretion.

Reuse Precaution Statement

Contents supplied STERILE (ethylene oxide (EO)). Do not use if sterile barrier is damaged. In the event of damaged packaging, call your representative of Taewoong Medical Co., Ltd. For a single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Storage: Store at room temperature (10~40°C).

Disposal Requirements: The introducer system of Uventa™ Urethral Stent must be properly sealed and disposed in compliance with the regulation of local or hospital at the end of its use.

***Warranty**

Taewoong Medical Co., Ltd. warrants that reasonable care has been used in the design and manufacture of this instrument.

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