# Niti-S & ComVi

## Enteral Colonic Stent

### Table of Contents

English.....	3
Français.....	8
Deutsch .....	13
Italiano.....	18
Español .....	23
Nederlands.....	28
Português .....	33
Türkçe.....	38
Русский .....	43
Čeština .....	48
Slovensky.....	53
Српски.....	58
Dansk.....	63
Ελληνικά.....	68
Lietuvių.....	73
Polski .....	78
Română.....	83
Svenska .....	88
Norsk .....	93
Български език .....	98
中國語.....	103

## ***MRI Information***

### **MR Conditional**

Non-clinical testing demonstrated that the Taewoong Medical Niti-S & ComVi Enteral Colonic Stent is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined, the Taewoong Medical Niti-S & ComVi Enteral Colonic Stent is expected to produce a maximum temperature rise of 2.2 °C after 15-minutes of continuous scanning.

### ***Artifact Information***

In non-clinical testing, the image artifact caused by the Taewoong Medical Niti-S & ComVi Enteral Colonic Stent extends approximately 5 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

## User's Manual

### 1. Description

The Niti-S & ComVi Enteral Colonic Stent consists of the implantable metallic stent and introducer system.

The stent is made of Nitinol wire. It is a flexible, fine mesh tubular prosthesis which has radiopaque markers on each end and at the center.

Model Name
Niti-S Enteral Colonic Uncovered Stent
Niti-S Enteral Colonic Covered Stent
ComVi Enteral Colonic Stent


**Figure1. Stent Model**

The Stent is loaded into an introducer system and upon deployment the stent imparts an outward radial force on the luminal surface of the colorectum to establish patency. There are two types of introducer; TTS and OTW. (See figure 2, 3)

Niti-S Full Covered Enteral Colonic Stents used in benign stricture can be removed; (see Warnings).


Full Covered Enteral Colonic stents can be repositioned after deployment (see Warnings).

#### A. TTS Type Introducer

**Figure 2. TTS Type Introducer System**

#### B. OTW Type Introducer

**Figure 3.OTW Type Introducer System  
(Distal release & Proximal release)**

TTS means Through The endoScope, and OTW means Over The Wire. The TTS type may be utilised for endoscopic delivery as the profile of the introducer system is compatible with a minimum 3.7mm working channel. The OTW type is not compatible with the endoscope due to the shorter length of the introducer system.

TTS type is recommended

- Where the target stricture region is narrow
- Where long-distance deploying is required
- Where exact and precise deployment is required

OTW is recommended

- Where fluoroscope is used

### 2. Principle of Operation

- For OTW Type; Distal release & TTS Type Introducer systems, The outer sheath is pulled back by immobilizing the hub in one hand, grasping the Y-connector with the other hand, and gently sliding the Y-connector along the pusher (2<sup>nd</sup> inner catheter) towards the hub. Retraction of the outer sheath releases the stent.
- For OTW Type; Proximal release systems, the distal part of outer sheath is pushed forward by immobilizing the Connector in one hand, grasping the hub with the other hand, and gently sliding the hub along the 2<sup>nd</sup> inner catheter towards the Connector. Retraction of the outer sheath releases the stent

### 3. Indication for Use

The Niti-S & ComVi Enteral Colonic Stent is intended for maintaining colon luminal patency in colon stricture caused by intrinsic and/or extrinsic malignant and/or benign stricture.

### WARRANTY

Taewoong Medical Co., LTD. warrants that reasonable care has been applied within the design and subsequent manufacturing process of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied

warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond Taewoong's control directly affect the instrument and the results obtained from its use. Taewoong's obligation under this warranty is limited to the repair or replacement of this instrument and Taewoong shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. Taewoong neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. Taewoong assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

#### **4. Contraindication**

The Niti-S & ComVi Enteral Colonic Stent is contraindicated for, but not limited to;

- Patient with ascites.
- Placement in polypoid lesions.
- Patient with bleeding disorder.
- Intra-abdominal abscess.
- Patients with coagulopathy.
- Strictures that do not allow passage of a guidewire.
- Any use other than those specifically outlined under indications for use.
- Removal or repositioning of fully deployed uncovered/bare Stents is contraindicated. (See Warnings).
- Suspected or impending perforation.
- Enteral ischemia.
- Recapturing a stent during its deployment is contraindicated.

#### **5. Warnings**

- The device should be used with caution and only after careful consideration in patients with elevated bleeding times, coagulopathies, or in patients with radiation colitis or proctitis.
- Chemoradiation therapy or radiotherapy alone may lead to tumor shrinkage and subsequent stent migration.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- Do not expose the introducer system to organic solvent (e.g. Alcohol)
- Do not use with Ethiodol or Lipiodol contrast media.
- Niti-S Full Covered Stents can be removed where the treating doctor's clinical assessment of the stented stricture indicates the Stent can be safely removed. Caution should be exercised in deciding to and when removing the Stent.
- Niti-S Full Covered Stent cannot be removed when there is tumor in-growth/over-growth/occlusion of the Stent lumen.
- Full Covered Stent can be repositioned immediately after deployment.
- Uncovered/bare Stents should not be removed once fully deployed, (See Contraindications).
- Do not attempt to recapture/reload a stent once its deployment is advanced.
- Niti-S Enteral Colonic Uncovered Stent and Niti-S Enteral Colonic Covered Stent (loaded in an Endoscopic Introducer) cannot be recaptured if the Y-connector has been pulled beyond the pusher's marker. Recapturing the stent in tortuous anatomy may damage the device. Recapturing more than twice may also cause damages to the silicone membrane and/or the stent wire.

#### **6. Potential complications**

Potential complications associated with the use of the colon stent may include, but are not limited to:

##### Procedural Complications

- Bleeding
- Stent misplacement or inadequate expansion
- Pain
- Death (other than that due to normal disease progression)
- Intestinal Perforation

##### Post Stent Placement Complications

- Bleeding
- Pain
- Perforation
- Bowel impaction
- Stent misplacement or migration
- Stent occlusion
- Tumor overgrowth
- Tumor ingrowth
- Stent fracture
- Fever
- Foreign body sensation
- Death (other than that due to normal disease progression)
- Sepsis
- Constipation
- Diarrhea

- Infection
- Peritonitis
- Ulceration
- Symptoms of tenesmus or urgency/incontinence

## 7. Equipment required

- Fluoroscope and/or Endoscope (with a channel size of 3.7 mm or greater)
- 0.035 in/ 0.89 mm guidewire (TTS)
- 0.038 in/ 0.97 mm guidewire (OTW; Distal & Proximal release)
- Introducer sheath appropriately sized for stent and introducer system

## 8. Precautions

Read the entire User's Manual thoroughly before using this device. It should only be used by or under the supervision of physicians thoroughly trained in the placement of stents. A thorough understanding of the techniques, principles, clinical applications and risks associated with this procedure is necessary before using the device.

- Care should be taken when removing the introducer system and guidewire immediately after stent deployment since this may result in stent dislodgement if the stent has not been adequately deployed.
- Care should be taken when performing dilation after the Stent has been deployed as this may result in perforation, bleeding, Stent dislodgement or Stent migration.
- The packaging and the device should be inspected prior to use.
- Use of fluoroscopy is recommended to ensure correct placement of the device.
- Check the expiration date "Use by". Do not use the device beyond the use by date.
- The Niti-S & ComVi Stent is supplied sterile. Do not use if the packaging is opened or damaged.
- The Niti-S & ComVi Stent is intended for single use only. Do not resterilize and/or reuse the device.

## 9. Instructions in the event of Damage

**WARNING:** Visually inspect the system for any sign of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this precaution may result in patient injury.

## 10. Procedure

### ① Examine stricture endoscopically and fluoroscopically

- Carefully examine both the proximal and distal segment of stricture endoscopically and/or fluoroscopically.
- The Internal luminal diameter should be measured exactly with endoscope and/or fluoroscope.

### ② Stent Size Determination

- Measure the length of the target stricture.
- Select a stent size that is 20 to 40mm longer than the measured length of the stricture in order to cover fully both ends of the lesion.
- Measure the diameter of the reference stricture - it is necessary to select a stent which has an unconstrained diameter about 1 to 4mm larger than the largest reference target diameter, to achieve secure placement.

### ③ Stent Deployment Preparation

- The Niti-S & ComVi Stent can be placed with the aid of fluoroscopy, and/or endoscopy.
- Pass a 0.035" (0.89 mm) guidewire to the level of the stricture. (TTS)
- Pass a 0.038" (0.97 mm) guidewire to the level of the stricture. (OTW)

#### A. TTS procedure

- Under the endoscopic guidance, insert an endoscope to the level of the obstruction, then introduce the guide wire through the working channel of the endoscopy. Advance until the guide wire across the target stricture to where the stent introducer system will be placed over the guide wire.
- Remove the stylet from the distal end of the introducer.
- Ensure that the valve of Y-connector connecting the inner sheath and outer sheath is locked by rotation proximal valve end in a clockwise direction to prevent premature stent deployment.
- Flush the inner lumen of introducer system.

#### B. OTW procedure

- Remove the stylet from the distal end of the introducer.
- Under the fluoroscopy guidance, insert a guide wire across the stricture to where the stent introducer system will be placed over the guide wire.
- Ensure that the valve of Y-connector connecting the inner sheath and outer sheath is locked by rotation proximal valve end in a clockwise direction to prevent premature stent deployment.

#### ④ Stent Deployment Procedure


Figure 4

**PRECAUTION:** Do not twist introducer system or employ a boring motion during the deployment as this may affect positioning and ultimate function of stent

##### A. OTW (Distal Release System) & TTS Introducer System

- Under the fluoroscope and/or endoscopic guidance, position the introducer system to the center of the target stricture exactly.
- Once the introducer system is in the correct position for deployment, unlock the proximal valve of the Y-connector by turning the valve more than twice in an anti-clockwise direction.
- To begin stent deployment, immobilize the hub in one hand and grasp the Y-connector with the other hand. Gently slide the Y-connector back along the pusher towards the hub.
- When the center X-ray marker reaches the center of target stricture, continue pulling back on the Y-connector until the stent is fully deployed. (See figure 4, 5)


Figure 5

**CAUTION** Do not push forward or pull backward on the hub with the stent partially deployed. The hub must be securely immobilized. Inadvertent movement of the hub may cause misalignment of the stent and possible damage to the target or stricture.

##### B. OTW (Proximal Release System)

- Under the fluoroscope and/or endoscopy guidance, position the introducer system to the center of target stricture exactly.
- Once the introducer system is in the correct position for deployment, unlock the proximal valve of the Connector by turning the valve more than twice in an anti-clockwise direction. The stent is now ready for deployment
- To begin stent deployment, immobilize the Connector in one hand and grasp the hub with the other hand. Gently slide the hub forward along the 2<sup>nd</sup> inner catheter towards the Connector.
- When the center X-ray marker reaches the center of target stricture, continue forwarding toward the connector until the stent is fully deployed. (See figure 4, 6)


Figure 6

**CAUTION** Do not push forward or pull backward on the Connector with the stent partially deployed. The Connector must be securely immobilized. Inadvertent movement of the Connector may cause misalignment of the stent and possible damage to the colorectum.

#### ⑤ After Stent Deployment

- Examine the stent fluoroscopically and/or endoscopically to confirm expansion.
- Carefully remove the introducer system, guidewire and endoscope from the patient. If excessive resistance is felt during removal, wait 3~5 minutes to allow further stent expansion (Place the inner sheath back into the outer sheath as the original state prior to removal.)
- Balloon dilatation inside the stent can be performed if on demand.

#### 11. Perform routine post implant procedure

- Assess the size and stricture of the Stent lumen. A Stent may require up to 1 to 3 days to expand fully.
- Doctor's experience and discretion can determine the appropriate drug regimen for each patient.
- After implantation, patient should remain on a soft diet until otherwise determined by the treating doctor.
- Observe the patient for development of any complications.

## **12. Instructions for removal of Niti-S Full Covered Stents (see Warnings)**

Visually examine the Stent for any tumor in-growth/over-growth into the Stent lumen or whether the Stent is occluded. If the Stent lumen is clear, carefully remove using a forcep and/or snare. Grasp the retrieval string and/or collapse the proximal end of the Stent then carefully retrieve the Stent. If the Stent cannot be easily withdrawn, do not remove the Stent.

**Caution:** Do not allow excessive force to remove the stent as it may cause disconnect to the retrieval string.

To reposition a Niti-S Full Covered Stent immediately after deployment, use forceps or a snare to grasp the retrieval string and gently adjust to the correct placement.

**Please note:** the stent can only be repositioned and/or removed proximally.

### **Reuse Precaution Statement**

Contents supplied STERILE (ethylene oxide (EO)). Do not use if sterile barrier is damaged. In the event of damaged packaging, call your Taewoong Medical Co., Ltd. representative. For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and /or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**Storage: Store at room temperature (10-40°C).**

**Disposal Requirements:** The introducer system of Niti-S & ComVi Enteral Colonic Stent must be properly sealed and disposed in compliance with the regulation of local or hospital at the end of its use.

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