Patient Information Leaflet for

PULSTA™ Transcatheter Pulmonary Valve System

Implant Description

The PULSTA heart valve is a medical device made from pig heart tissue, and has a stent frame (nickel-titanium). It helps replace a damaged heart valve that can no longer properly control blood flow to the lungs. Doctors can insert it using a thin tube (catheter) without open-heart surgery.

This heart valve is designed for those who have any of the following conditions.

- Narrowing of the pulmonary valve
- Leaking of the pulmonary valve
- Problems with blood flow from the right side of the heart to the lungs

Information on safe use

- Patient Implant Card

After the procedure, you will get an implant card with information about your PULSTA heart valve. Carry this card with you at all times and show it before any dental treatment (such as teeth scaling) or medical procedure. It's important to tell your doctors about your heart valve. The doctor may need to give you antibiotics before a procedure.

- Magnetic Resonance Imaging (MRI) scan

Before getting an MRI, always show your implant card to your doctor to make sure the scan is safe for your heart valve. The heart valve can be safely scanned under certain MRI conditions. For more details, refer to the IFU at www.taewoongmedical.com/manual, or send an e-mail to contact@stent.net.

- Metal Detector

Because the implanted valve might set off metal detectors, always carry your implant card and tell airport security.

- Allergic Reactions

Allergic reactions to metal (nickel) or valve materials may occur. If you notice symptoms such as rash, swelling, or difficulty breathing, call your doctor right away.

- Infection Risks

People with infections (like endocarditis) should not use this device. If you have symptoms of infection such as fever or extreme tiredness, call your doctor right away.

Follow-up Care

- After The Procedure

After the procedure, doctors will check your heart valve in the recovery room and monitor for any issues. Most patients go home within several days, and light walking is usually fine 24 hours after the procedure.

However, the recovery time varies for each person, so it's important to follow your doctor's advice.

You may need to take blood thinners or other medications to prevent blood clots. Be sure to take all medications exactly as prescribed. Do not change or stop your medication without consulting your doctor first.

- After Discharge

Call your doctor right away if you have bleeding, pain, trouble breathing, fever, extreme tiredness, or an irregular heartbeat. These symptoms could be a sign of a serious problem. It is very important to get medical care quickly if an infection or other issues are suspected.

You must visit the hospital regularly to check your heart valve and overall health. Heart ultrasound tests are needed, and your doctor may also perform extra tests such as CT or MRI scans, if necessary.

Device Lifetime

The heart valve was tested in a lab for 5 years without any problems, but its expected lifetime varies for each patient, and it may need to be replaced by a doctor's decision when the heart valve shows significant dysfunction.

Material Information

The PULSTA™ TPV is made of the following materials, and contacts the patient.

Section		Material and Substance			Mass ¹⁾
Valve leaflets		Porcine pericardium			1.0884 g
Stent frame	Nitinol stent	Nickel-titanium alloy	Nickel	7440-02-0	1.0206 g
			Carbon	7440-44-0	
			Cobalt	7440-48-4	
			Copper	7440-50-8	
			Chromium	7440-47-3	
			Hydrogen	1333-74-0	
			Iron	7439-89-6	
			Niobium	7440-03-1	
			Nitrogen	7727-37-9	
			Oxygen	7782-44-7	
			Titanium	7440-32-6	
	Radiopaque marker	Platinum-iridium alloy	Platinum	7440-06-4	0.0035 g
			Iridium	7439-88-5	
	Connecting tube	Stainless steel 316L ²⁾	Iron	7439-89-6	0.0045 g
			Carbon	7440-44-0	

Section		Material and Substance			Mass ¹⁾
			Manganese	7439-96-5	
			Phosphorus	7723-14-0	
			Sulfur	7704-34-9	
			Silicon	7440-21-3	
			Chromium	7440-47-3	
			Nickel	7440-02-0	
			Molybdenum	7439-98-7	
			Nitrogen	7727-37-9	
			Cobalt	7440-48-4	
Polyester suture		Polyethylene terephthalate			0.0645 g

- 1) The mass of each material was measured based on the largest valve size (TPV3238).
- 2) This stainless steel 316L contains the following substance classified as CMR 1B in a concentration above 0.1% weight by weight:
Cobalt (CAS No. 7440-48-4; EC No. 231-158-0) Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not pose an increased risk of cancer or adverse reproductive effects.

Implant Card Symbols

	Patient name or patient ID
	Date of implantation
	Name and address of the implanting healthcare institution/provider
	Information website for patients
	Device name
	Unique device identifier
	Model number
	Serial number
	Manufacturer
	MR conditional: An MRI scan is safe only under specific conditions. Consult your doctor

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Taewoong Medical Co., Ltd.
Emergo Consulting (UK) Limited