




Through The Scope (TTS) Esophageal Stent



“The practical solution for tight, narrow or tortuous anatomies with 10.5 Fr delivery system”

MANUFACTURER : TAEWOONG MEDICAL Co., Ltd 14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, South Korea 10022
TEL: +82 31 996 0641~4 FAX: +82 31 996 0645~6

SALES OFFICE 4801 Wilshire Blvd. Suite 280 Los Angeles, CA 90010
TEL: 323 692 0028 TOLL FREE: 844 303 1112 FAX: 323 692 0038 E-MAIL: info@taewoongusa.com

U.S. Distributor: Cook Medical Inc.



Through The Scope (TTS) Esophageal Stent

● Easy and Simple stenting through the scope channel

Preloaded in a 10.5Fr delivery system for esophageal fully covered/partially covered stent


The practical solution for tight, narrow or tortuous anatomies

Retrieval string may help removal of the stent during the initial placement procedure

Silicone coating designed to resist tissue in-growth

Clearer yellow marker is designed to be accurate placement under endoscopic visualization

Radiopaque marker : 4(four) at both ends & 2(two) in the middle



● Ordering Information

TTS (Through The Scope)

► Fully covered

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)
EST1806F	18	6	10.5	160 180 (USA)
EST1808F		8		
EST1810F		10		
EST1812F		12		
EST1814F		14		
EST1815F		15		
EST2006F	20	6	10.5	160 180 (USA)
EST2008F		8		
EST2010F		10		
EST2012F		12		
EST2014F		14		
EST2015F		15		
EST1806F-18	18	6	10.5	180
EST1808F-18		8		
EST1810F-18		10		
EST1812F-18		12		
EST1814F-18		14		
EST1815F-18		15		
EST2006F-18	20	6	10.5	180
EST2008F-18		8		
EST2010F-18		10		
EST2012F-18		12		
EST2014F-18		14		
EST2015F-18		15		

► Partially covered

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)
EST1806B	18	6	10.5	160 180 (USA)
EST1808B		8		
EST1810B		10		
EST1812B		12		
EST1814B		14		
EST1815B		15		
EST2006B	20	6	10.5	160 180 (USA)
EST2008B		8		
EST2010B		10		
EST2012B		12		
EST2014B		14		
EST2015B		15		
EST1806B-18	18	6	10.5	180
EST1808B-18		8		
EST1810B-18		10		
EST1812B-18		12		
EST1814B-18		14		
EST1815B-18		15		
EST2006B-18	20	6	10.5	180
EST2008B-18		8		
EST2010B-18		10		
EST2012B-18		12		
EST2014B-18		14		
EST2015B-18		15		

* Available in USA (The sizes marked in blue)

* EST**_F-22 / EST**_B-22: Available Delivery system length 220cm

