




S Esophageal Stent



MANUFACTURER : TAEWOONG MEDICAL Co., Ltd 14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, South Korea 10022
TEL: +82 31 996 0641~4 FAX: +82 31 996 0645~6

SALES OFFICE 4801 Wilshire Blvd. Suite 280 Los Angeles, CA 90010
TEL: 323 692 0028 TOLL FREE: 844 303 1112 FAX: 323 692 0038 E-MAIL: info@taewoongusa.com

U.S. Distributor: Cook Medical Inc.



Niti-S™ S Esophageal Stent

for malignant esophageal strictures

● Features

- Fixed cell with braided construction
 - High flexibility and optimal radial force
 - Both head ends (8mm larger than trunk) help to minimize migration
- Silicone covering and soft round ends
 - Reduce tissue ingrowth and hyperplasia reaction
- Visible green suture for easy removal
- Radiopaque marker : 4(four) at both ends & 2(two) in the middle



● Ordering Information

OTW (Over The Wire)

► Fully covered

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)
ES1606F	16	6	16	70
ES1608F		8		
ES1610F		10		
ES1612F		12		
ES1614F		14		
ES1615F		15		
ES1806F	18	6	16	
ES1808F		8		
ES1810F		10		
ES1812F		12		
ES1814F		14		
ES1815F		15		
ES2006F	20	6	20	
ES2008F		8		
ES2010F		10		
ES2012F		12		
ES2014F		14		
ES2015F		15		

► Partially covered

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)
ES1606B	16	6	16	70
ES1608B		8		
ES1610B		10		
ES1612B		12		
ES1614B		14		
ES1615B		15		
ES1806B	18	6	16	
ES1808B		8		
ES1810B		10		
ES1812B		12		
ES1814B		14		
ES1815B		15		
ES2006B	20	6	20	
ES2008B		8		
ES2010B		10		
ES2012B		12		
ES2014B		14		
ES2015B		15		

Released Article

* Fully covered, retrievable self-expanding metal stents(Niti-s) in palliation of malignant dysphagia: Long-term results of a prospective study
by Sung Jun Choi et al [Scand J Gastroenterol. 2011 Jul;46(7-8):875-80]