

ORDERING INFORMATION

Distal Release					Proximal Release								
Code	Stent		Delivery		Code	Stent		Delivery					
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)		Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)				
ES1606FD	16	6	18	70	ESP1606FD	16	6	18	70				
ES1608FD		8			ESP1608FD		8						
ES1610FD		10			ESP1610FD		10						
ES1612FD		12			ESP1612FD		12						
ES1614FD		14			ESP1614FD		14						
ES1615FD	15	ESP1615FD			15								
ES1806FD	18	6			20	ESP1806FD	18			6	20		
ES1808FD		8				ESP1808FD				8			
ES1810FD		10				ESP1810FD				10			
ES1812FD		12				ESP1812FD				12			
ES1814FD		14				ESP1814FD				14			
ES1815FD	15	ESP1815FD				15							
ES2006FD	20	6				22	ESP2006FD			20		6	22
ES2008FD		8					ESP2008FD					8	
ES2010FD		10					ESP2010FD					10	
ES2012FD		12	ESP2012FD				12						
ES2014FD		14	ESP2014FD				14						
ES2015FD	15	ESP2015FD	15										
ES2206FD	22	6	24				ESP2206FD	22		6		24	
ES2208FD		8					ESP2208FD			8			
ES2210FD		10					ESP2210FD			10			
ES2212FD		12			ESP2212FD		12						
ES2214FD		14			ESP2214FD		14						
ES2215FD	15	ESP2215FD			15								
ES2406FD	24	6			28		ESP2406FD	24		6	28		
ES2408FD		8					ESP2408FD			8			
ES2410FD		10					ESP2410FD			10			
ES2412FD		12				ESP2412FD	12						
ES2414FD		14				ESP2414FD	14						
ES2415FD	15	ESP2415FD				15							
ES2806FD	28	6		28		ESP2806FD	28	6	28				
ES2808FD		8				ESP2808FD		8					
ES2810FD		10				ESP2810FD		10					
ES2812FD		12	ESP2812FD			12							
ES2814FD		14	ESP2814FD			14							
ES2815FD	15	ESP2815FD	15										

*Stents with retrieval strings on both ends: ES**##FDR2

*Stents with retrieval strings on both ends: ESP**##FDR2

RELEASED ARTICLE

- * Double-layered covered stent for the treatment of malignant oesophageal obstructions: Systematic review and meta-analysis by Zeiad Hussain et al [World J Gastroenterol. 2016 Sep 14;22(34):7841-50]
- * Double-type metallic stents efficacy for the management of post-operative fistulas, leakages, and perforations of the upper gastrointestinal tract by Jean-Michel Gonzalez et al [Surg Endosc 2015;29:2013-2018]
- * Outcomes Following Oesophageal Stent Insertion for Palliation of Malignant Strictures: A Large Single Center Series by Nicholas J. Battersby, BMEDSc, MARCS et al [J Surg Oncol. 2012 Jan;105(1):60-5]
- * New Design Esophageal Stents for the Palliation of Dysphagia From Esophageal or Gastric Cardia Cancer: A Randomized Trial by Els M.L. Verschuur, M.Sc et al [Am J Gastroenterol. 2008 Feb;103(2):304-12]
- * A new esophageal stent design (Niti-S stent) for the prevention of migration: a prospective study in 42 patients by NICHOLAS J. BATTERSBY, BMedSc, MRCS et al [Gastrointest Endosc. 2006 Jan;63(1):134-40]