


CONIO™ II Esophageal Stent



**Optimal Solution for Benign and Malignant
Upper, Middle Esophageal Strictures**

HEAD OFFICE

14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, Korea 10022

SALES OFFICE

Bldg.C 4F, 158 Haneulmaeul-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do, Korea 10355
E: contact@stent.net T: +82 31 904 6153 F: +82 31 904 6157

www.taewoongmedical.com

Niti-S™ CONIO™ II Esophageal Stent

for benign and malignant upper, middle esophageal stricture

Features



HEAD

- Silicone-covered S type for better palliation
- 2mm bigger than the body

BODY

- Conformable D type with PTFE membrane to prevent stent-related complications

BUMPER

- Silicone-coated Double Ring for anti-migration
- 6mm bigger than the body at the peak
- The diameter of the hole is 2mm

Ordering Information

Distal Release				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)
EK1008FN	10	8	14	70
EK1010FN		10		
EK1012FN		12		
EK1208FN	12	8	16	
EK1210FN		10		
EK1212FN		12		
EK1408FN	14	8	18	
EK1410FN		10		
EK1412FN		12		
EK1608FN	16	8	18	
EK1610FN		10		
EK1612FN		12		