

ORDERING INFORMATION

Distal Release

► A retrieval string on the proximal end

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ES2218F	22	18	20	70
ES2223F		23		
ES2418F	24	18		
ES2423F		23		
ES2818F	28	18		
ES2823F		23		

► Retrieval string on both ends

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ES2218FR2	22	18	20	70
ES2223FR2		23		
ES2418FR2	24	18		
ES2423FR2		23		
ES2818FR2	28	18		
ES2823FR2		23		

RELEASED ARTICLE

- * A novel management of post-oesophagectomy gastro-pleural fistula
by Javaid Ishtiaq et al. [J gastrointest Oncol 2016;7(6):E93-E97]
- * An Endoscopic Strategy Combining Mega Stents and Over-The-Scope Clips for the Management of Post-Bariatric Surgery Leaks and Fistula by Hany M. Shehab et al [Obes Surg. 2016 May;26(5):941-8]
- * Staple-line leak after sleeve gastrectomy in obese patients: A hot topic in bariatric surgery
by Giuseppe Galloro et al [World J Gastrointest Endosc. 2015 Jul 25;7(9):843-6]
- * Use of sleeve-customized self-expandable metal stents for the treatment of staple-line leakage after laparoscopic sleeve gastrectomy by Sigal Fishman, MD et al [Gastrointest Endosc. 2015 May;81(5):1291-4]
- * A novel dedicated endoscopic stent for staple line leaks after laparoscopic sleeve gastrectomy: a case series
by Giuseppe Galloro, MD et al [Surg Obes Relat Dis. 2014 Jul-Aug;10(4):607-11]
- * Mega stent: a new option for management of leaks following laparoscopic sleeve gastrectomy
by Jahangeer Basha et al [Endoscopy. 2014;46 Suppl 1 UCTN:E49-50]