


Optimos™ Polypectomy Snare

360° Rotation

The rotatable locking controller has been designed to provide a 360° rotation of the snare loop which allows accurate access to polyps from any angle.



HEAD OFFICE

14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, Korea 10022

SALES OFFICE

Bldg.C 4F, 158 Haneulmaeul-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do, Korea 10355
E: contact@stent.net T: +82 31 904 6153 F: +82 31 904 6157

www.taewoongmedical.com

TaeWoong
MEDICAL 

Ordering Information

360° Rotatable

Art. No.	Type	Loop	Outer Sheath		Wire				
		Diameter (mm)	Diameter (mm)	Working Length (mm)	Diameter (mm)				
KSN16R-10	Oval	10	2.4	1650	0.45				
KSN23R-10				2300					
KSN16R-15				15		1650			
KSN23R-15						2300			
KSN16R-20						20	1650		
KSN23R-20							2300		
KSN16R-25		25					1650		
KSN23R-25							2300		
KSN16R-30				30			1650		
KSN23R-30							2300		
KSN23C-10						Round Type	2.4	2300	0.45
KSN23C-15									
KSN23C-20	16								
KSN23C-25	20								
KSN23C-30	25								
KSN23C-30	30								

Non Rotatable

Art. No.	Type	Loop	Outer Sheath		Wire				
		Diameter (mm)	Diameter (mm)	Working Length (mm)	Diameter (mm)				
KSN16RN-10	Oval	10	2.4	1650	0.45				
KSN23RN-10				2300					
KSN16RN-15				15		1650			
KSN23RN-15						2300			
KSN16RN-20						20	1650		
KSN23RN-20							2300		
KSN16RN-25		25					1650		
KSN23RN-25							2300		
KSN16RN-30				30			1650		
KSN23RN-30							2300		
KSN23RNS-10						Oval	2.4	2300	0.3
KSN23RNS-15									
KSN23RNS-20	15								
KSN23RNS-20	20								

Optimos™
Polypectomy Snare



360° Rotation

The rotatable locking controller has been designed to provide a 360° rotation of the snare loop which allows accurate access to polyps from any angle.