


# Optimos™ Injector



21G



23G



25G

## Optimos™ Injector

- Ergonomic handle design gives stability and helps to have a precise needle control
- Locking system prevents the handle to get loose
- Wider diameter of the catheter facilitates drugs injection
- Different colored handles help users to identify the appropriate needle size

### HEAD OFFICE

14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, Korea 10022

### SALES OFFICE

Bldg.C 4F, 158 Haneulmaeul-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do, Korea 10355  
E: [contact@stent.net](mailto:contact@stent.net) T: +82 31 904 6153 F: +82 31 904 6157

[www.taewoongmedical.com](http://www.taewoongmedical.com)

**TaeWoong**  
MEDICAL 



## Features

- Ergonomic handle design gives stability and helps to have a precise needle control
- Locking system prevents the handle to get loose
- Wider diameter of the catheter facilitates drugs injection
- Different colored handles help users to identify the appropriate needle size

## Ordering Information

Item no.	Needle		Sheath		Handle color
	diameter	length	diameter	length	
OIP1604-21G	21G	4mm	2.4mm	1650mm	Light Blue
OIP1604-23G	23G	4mm	2.4mm	1650mm	Purple
OIP1604-25G	25G	4mm	2.4mm	1650mm	Yellow
OIP2304-21G	21G	4mm	2.4mm	2300mm	Light Blue
OIP2304-23G	23G	4mm	2.4mm	2300mm	Purple
OIP2304-25G	25G	4mm	2.4mm	2300mm	Yellow

# Optimos™ Injector

## Different colored handles

*help users to identify the appropriate needle size*

