




Optimos™

STONE BASKET



Optimos™ STONE BASKET

- **Optimos™ Stone basket** is for endoscopic removal of biliary/pancreatic stones and foreign bodies
- The basket designed to provide rotation allow users capture stones well
- The spiral type with eight wires allow users capture and retain stones for an easy removal
- The separated handle offers good environment for mechanical lithotripsy in case of emergency



FOLD TYPE



SPIRAL TYPE

HEAD OFFICE

14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, Korea 10022


SALES OFFICE

Bldg.C 4F, 158 Haneulmaeul-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do, Korea 10355
E: contact@stent.net T: +82 31 904 6153 F: +82 31 904 6157

www.taewoongmedical.com

Art. No.	Type	Basket	Sheath		Wire
		Opening Width (mm)	Diameter (mm)	Working Length (mm)	
SB1825F4-15	Fold Angular	15	1.8	2500	4 wires
SB1825F4-20		20			
SB2425F4-20		20	2.4	2500	4 wires
SB2425F4-30		30			
SB2425F4-40		40			
SB2425F8-20		20			
SB2425F8-30		30			
SB1821T8-20	Spirality	20	1.8	2100	8 wires
SB1821T8-30		30			
SB2425T8-20		20	2.4	2500	8 wires
SB2425T8-30		30			

O P T I M O S S T O N E B A S K E T



- Optimos™ Stone basket is for endoscopic removal of biliary/pancreatic stones and foreign bodies.
- The basket designed to provide a rotation allow users capture stones well
- The spiral type with eight wires allow users capture and retain stones for an easy removal
- The separated handle offers good environment for mechanical lithotripsy in case of emergency