



Ordering information



Pulsta™ Transcatheter Pulmonary Valve System

| Model name | Outer diameter | Flare diameter | Total length | Applicable delivery system model name |
|------------|----------------|----------------|--------------|---------------------------------------|
| TPV1828 | Ø18 mm | Ø22 mm | 28 mm | DSF1828 |
| TPV1838 | | | 38mm | DSF1838 |
| TPV2028 | Ø20 mm | Ø24 mm | 28 mm | DSF2028 |
| TPV2038 | | | 38 mm | DSF2038 |
| TPV2231 | Ø22 mm | Ø26 mm | 31 mm | DSF2231 |
| TPV2238 | | | 38 mm | DSF2238 |
| TPV2431 | Ø24 mm | Ø28 mm | 31 mm | DSF2431 |
| TPV2438 | | | 38 mm | DSF2438 |
| TPV2633 | Ø26 mm | Ø30 mm | 33 mm | DSF2633 |
| TPV2638 | | | 38 mm | DSF2638 |
| TPV2833 | Ø28 mm | Ø32 mm | 33 mm | DSF2833 |
| TPV2838 | | | 38 mm | DSF2838 |
| TPV3038 | Ø30 mm | Ø34 mm | 38 mm | DSF3038 |
| TPV3238 | Ø32 mm | Ø36 mm | 38 mm | DSF3238 |



Delivery system

| Model name | Outer diameter | Usable length | Total length |
|------------|----------------|---------------|--------------|
| DSF1828 | 18 Fr | 110 cm | 138 cm |
| DSF1838 | | | |
| DSF2028 | | | |
| DSF2038 | | | |
| DSF2231 | | | |
| DSF2238 | | | |
| DSF2431 | | | |
| DSF2438 | | | |
| DSF2633 | | | |
| DSF2638 | | | |
| DSF2833 | | | |
| DSF2838 | 20 Fr | | |
| DSF3038 | | | |
| DSF3238 | | | |