


Through The Scope (TTS) Delivery System

S Esophageal Stent (Covered)

“The practical solution for tight, narrow or tortuous anatomies with 10.5 Fr delivery system”



TTS Esophageal stent passing through the tight and narrow stricture

Dr. Giorgio Battaglia, Italy

HEAD OFFICE

14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, Korea 10022

SALES OFFICE

Bldg.C 4F, 158 Haneulmaeul-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do, Korea 10355
E: contact@stent.net T: +82 31 904 6153 F: +82 31 904 6157


www.taewoongmedical.com

Niti-S™ Through The Scope (TTS) Delivery System

S Esophageal Stent (Covered)

Easy and Simple stenting through the scope channel

- Preloaded in a 10.5Fr delivery system for esophageal fully covered/partially covered stent
- The practical solution for tight, narrow or tortuous anatomies
- Proximal retrieval string may help removal of the stent
- Silicone coating designed to resist tissue in-growth
- Yellow marker on the black inner sheath for accurate placement under endoscopic visualization
- Radiopaque marker: 4 (four) at both ends & 2 (two) in the middle



Ordering Information

TTS (Through The Scope)

► Full covered-type

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)
EST1806F	18	6	10.5	160
EST1808F		8		
EST1810F		10		
EST1812F		12		
EST1814F		14		
EST1815F		15		
EST2006F	20	6		
EST2008F		8		
EST2010F		10		
EST2012F		12		
EST2014F		14		
EST2015F		15		
EST2206F	22	6		
EST2208F		8		
EST2210F		10		
EST2212F		12		
EST2214F		14		
EST2215F		15		

► Both bare-type

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)
EST1806B	18	6	10.5	160
EST1808B		8		
EST1810B		10		
EST1812B		12		
EST1814B		14		
EST1815B		15		
EST2006B	20	6		
EST2008B		8		
EST2010B		10		
EST2012B		12		
EST2014B		14		
EST2015B		15		
EST2206B	22	6		
EST2208B		8		
EST2210B		10		
EST2212B		12		
EST2214B		14		
EST2215B		15		

*Stents on the 180cm delivery system: EST****F-18

*Stents on the 220cm delivery system: EST****F-22

*Stents with removal strings on both ends available only with the 160cm delivery system: EST****FR2

*Stents on the 180cm delivery system: EST****B-18

*Stents on the 220cm delivery system: EST****B-22