


DOUBLE™

Esophageal Stent (Anti-reflux)



HEAD OFFICE

14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, Korea 10022

SALES OFFICE

Bldg.C 4F, 158 Haneulmaeul-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do, Korea 10355
E: contact@stent.net T: +82 31 904 6153 F: +82 31 904 6157


www.taewoongmedical.com

Niti-S™ **DOUBLE™** Esophageal Stent (Anti-reflux)

for malignant esophageal stricture and Preventing gastroesophageal reflux

Features

- Double layered design with Anti-reflux skirt
 - PTFE skirt blocks gastric reflux with the stent placement at EG junction
 - Additional uncovered outer mesh helps to resist migration
- Retrieval string at proximal end helps repositioning
- Radiopaque marker: 4 (four) at both ends & 2 (two) in the middle



Ordering Information

Distal Release									
Code	Stent		Delivery		Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)		Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
EA1606FD	16	6	18	70	EA2206FD	22	6	20	70
EA1608FD		8			EA2208FD		8		
EA1610FD		10			EA2210FD		10		
EA1612FD		12			EA2212FD		12		
EA1614FD		14			EA2214FD		14		
EA1615FD		15			EA2215FD		15		
EA1806FD	18	6	20		EA2406FD	24	6	22	
EA1808FD		8			EA2408FD		8		
EA1810FD		10			EA2410FD		10		
EA1812FD		12			EA2412FD		12		
EA1814FD		14			EA2414FD		14		
EA1815FD		15			EA2415FD		15		
EA2006FD	20	6	20		EA2806FD	28	6	22	
EA2008FD		8			EA2808FD		8		
EA2010FD		10			EA2810FD		10		
EA2012FD		12			EA2812FD		12		
EA2014FD		14			EA2814FD		14		
EA2015FD		15			EA2815FD		15		