




# MEGA™ Esophageal Stent



**HEAD OFFICE**

14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, Korea 10022

**SALES OFFICE**

Bldg.C 4F, 158 Haneulmaeul-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do, Korea 10355  
E: [contact@stent.net](mailto:contact@stent.net) T: +82 31 904 6153 F: +82 31 904 6157

[www.taewoongmedical.com](http://www.taewoongmedical.com)




# Niti-S™ MEGA™ Esophageal Stent

for benign and malignant esophageal stricture

## Features

- Specially designed soft and flexible body
- Large diameter and long length of the stent: Prevent migration
- Full silicone covering allows easy removal
- Radiopaque marker: 4 (four) at both ends & 2 (two) in the middle



## Ordering Information

### Distal Release

► A retrieval string on the proximal end

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ES2218F	22	18	20	70
ES2223F		23		
ES2418F	24	18		
ES2423F		23		
ES2818F	28	18		
ES2823F		23		

► Retrieval string on both ends

Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ES2218FR2	22	18	20	70
ES2223FR2		23		
ES2418FR2	24	18		
ES2423FR2		23		
ES2818FR2	28	18		
ES2823FR2		23		