

Biliary Covered Stent [Bumpy]



HEAD OFFICE

14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, Korea 10022

SALES OFFICE

Bldg.C 4F, 158 Haneulmaeul-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do, Korea 10355
E: contact@stent.net T: +82 31 904 6153 F: +82 31 904 6157

www.taewoongmedical.com

Niti-S™ Biliary Covered Stent [Bumpy]

Intended for maintaining biliary luminal patency in malignant and/or benign strictures and benign pancreatic ductal strictures

Features

- Both flare ends and the irregular body cell size of segmental radial force prevent the risk of migration
- The short length of the stent reduces stent-related complications
- Platinum radiopaque marked, 10cm long retrieval string helps easy removal
- Radiopaque marker: 3 (three) at both ends & 2 (two) in the middle



Ordering Information

Endoscopic Approach				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (fr)	Usable Length (cm)
BK0604CW2	6	4	8.5	180
BK0606CW2		6		
BK0608CW2		8		
BK0804CW2	8	4		
BK0806CW2		6		
BK0808CW2		8		
BK1004CW2	10	4		
BK1006CW2		6		
BK1008CW2		8		