
Niti-S Biliary Covered Stent [Giobor]

for malignant biliary strictures



Features

- Partially covered design
- (1) The long covered portion (70%)
- (2) The bare portion (30%)
- Anti-migration flared end
- Radiopaque markers: Three (3) at both ends & Two (2) at the end of the covered part



Ordering Information

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BS0806FP	8	6	8.5	180
BS0808FP		8		
BS0810FP		10		
BS0812FP		12		
BS1006FP	10	6	8.5	180
BS1008FP		8		
BS1010FP		10		
BS1012FP		12		