


# Niti-S Biliary Uncovered Stent [T-Type / Y-Type]

for malignant biliary strictures



## Features

- **Exclusive Design**
- Larger central interstices for additional access and/or reintervention

# Ordering Information

Percutaneous				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TT1004A	10	4	8	50
TT1005A		5		
TT1006A		6		
TT1007A		7		
TT1008A		8		
TT1009A		9		
TT1010A		10		
TT1012A		12		

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BY1004A	10	4	8	180
BY1005A		5		
BY1006A		6		
BY1007A		7		
BY1008A		8		
BY1009A		9		
BY1010A		10		
BY1012A		12		