


Niti-S Esophageal Covered Stent [Full Covered-Type]

for malignant and/or benign stricture and tracheoesophageal fistula



Features

- **Conformable Stent Structure:** enables stent to conform to the esophagus movement
- **Radiopaque marker:** 4 (four) at both ends of the flanges & 2 (two) at both ends of the body of the stent

Ordering Information

Distal Release				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
EK1806F	18	8	16	70
EK1808F		10		
EK1810F		12		
EK1812F		14		
EK1814F		15		
EK1815F		16		
EK1818F		18		
EK1820F		20		
EK2006F	20	8	18	70
EK2008F		10		
EK2010F		12		
EK2012F		14		
EK2014F		15		
EK2015F		16		
EK2018F		18		
EK2020F		20		
EK2206F	22	8	18	70
EK2208F		10		
EK2210F		12		
EK2212F		14		
EK2214F		15		
EK2215F		16		
EK2218F		18		
EK2406F	24	6	18	70
EK2408F		8		
EK2410F		10		
EK2412F		12		
EK2414F		14		
EK2415F		15		
EK2806F	28	6	18	70
EK2808F		8		
EK2810F		10		
EK2812F		12		
EK2814F		14		
EK2815F		15		