


Niti-S Esophageal Covered Stent [Conio]

for malignant and/or benign stricture and tracheoesophageal fistula



Features

- **Silicone covering:** Reduce the risk of tumor ingrowth
- **Visible green suture for easy removal**
- **Radiopaque marker:** 4 (four) at both ends & 2 (two) in the middle

Ordering Information

Distal Release				
Distal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ES1006FN	10	6	12	70
ES1008FN		8		
ES1010FN		10		
ES1012FN		12		
ES1014FN		14		
ES1015FN		15		
ES1206FN	12	6	14	70
ES1208FN		8		
ES1210FN		10		
ES1212FN		12		
ES1214FN		14		
ES1215FN		15		
ES1406FN	14	6	14	70
ES1408FN		8		
ES1410FN		10		
ES1412FN		12		
ES1414FN		14		
ES1415FN		15		
ES1606FN	16	6	16	70
ES1608FN		8		
ES1610FN		10		
ES1612FN		12		
ES1614FN		14		
ES1615FN		15		
ES1806FN	18	6	16	70
ES1808FN		8		
ES1810FN		10		
ES1812FN		12		
ES1814FN		14		
ES1815FN		15		

Distal Release				
Distal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ES2006FN	20	6	20	70
ES2008FN		8		
ES2010FN		10		
ES2012FN		12		
ES2014FN		14		
ES2015FN		15		
ES2206FN	22	6	20	70
ES2208FN		8		
ES2210FN		10		
ES2212FN		12		
ES2214FN		14		
ES2215FN		15		
ES2406FN	24	6	22	70
ES2408FN		8		
ES2410FN		10		
ES2412FN		12		
ES2414FN		14		
ES2415FN		15		
ES2806FN	28	6	22	70
ES2808FN		8		
ES2810FN		10		
ES2812FN		12		
ES2814FN		14		
ES2815FN		15		

Ordering Information

Proximal Release				
Distal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ESP1006FN	10	6	14	70
ESP1008FN		8		
ESP1010FN		10		
ESP1012FN		12		
ESP1014FN		14		
ESP1015FN		15		
ESP1206FN	12	6	14	70
ESP1208FN		8		
ESP1210FN		10		
ESP1212FN		12		
ESP1214FN		14		
ESP1215FN		15		
ESP1406FN	14	6	14	70
ESP1408FN		8		
ESP1410FN		10		
ESP1412FN		12		
ESP1414FN		14		
ESP1415FN		15		
ESP1606FN	16	6	16	70
ESP1608FN		8		
ESP1610FN		10		
ESP1612FN		12		
ESP1614FN		14		
ESP1615FN		15		
ESP1806FN	18	6	16	70
ESP1808FN		8		
ESP1810FN		10		
ESP1812FN		12		
ESP1814FN		14		
ESP1815FN		15		

Proximal Release				
Distal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ESP2006FN	20	6	20	70
ESP2008FN		8		
ESP2010FN		10		
ESP2012FN		12		
ESP2014FN		14		
ESP2015FN		15		
ESP2206FN	22	6	20	70
ESP2208FN		8		
ESP2210FN		10		
ESP2212FN		12		
ESP2214FN		14		
ESP2215FN		15		
ESP2406FN	24	6	22	70
ESP2408FN		8		
ESP2410FN		10		
ESP2412FN		12		
ESP2414FN		14		
ESP2415FN		15		
ESP2806FN	28	6	22	70
ESP2808FN		8		
ESP2810FN		10		
ESP2812FN		12		
ESP2814FN		14		
ESP2815FN		15		

Ordering Information

Distal Release				
Distal & Proximal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ES1006FNR2	10	6	12	70
ES1008FNR2		8		
ES1010FNR2		10		
ES1012FNR2		12		
ES1014FNR2		14		
ES1015FNR2		15		
ES1206FNR2	12	6	14	70
ES1208FNR2		8		
ES1210FNR2		10		
ES1212FNR2		12		
ES1214FNR2		14		
ES1215FNR2		15		
ES1406FNR2	14	6	14	70
ES1408FNR2		8		
ES1410FNR2		10		
ES1412FNR2		12		
ES1414FNR2		14		
ES1415FNR2		15		
ES1606FNR2	16	6	16	70
ES1608FNR2		8		
ES1610FNR2		10		
ES1612FNR2		12		
ES1614FNR2		14		
ES1615FNR2		15		
ES1806FNR2	18	6	16	70
ES1808FNR2		8		
ES1810FNR2		10		
ES1812FNR2		12		
ES1814FNR2		14		
ES1815FNR2		15		

Distal Release				
Distal & Proximal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ES2006FNR2	20	6	20	70
ES2008FNR2		8		
ES2010FNR2		10		
ES2012FNR2		12		
ES2014FNR2		14		
ES2015FNR2		15		
ES2206FNR2	22	6	20	70
ES2208FNR2		8		
ES2210FNR2		10		
ES2212FNR2		12		
ES2214FNR2		14		
ES2215FNR2		15		
ES2406FNR2	24	6	22	70
ES2408FNR2		8		
ES2410FNR2		10		
ES2412FNR2		12		
ES2414FNR2		14		
ES2415FNR2		15		
ES2806FNR2	28	6	22	70
ES2808FNR2		8		
ES2810FNR2		10		
ES2812FNR2		12		
ES2814FNR2		14		
ES2815FNR2		15		

Ordering Information

Proximal Release				
Distal & Proximal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ESP1006FNR2	10	6	14	70
ESP1008FNR2		8		
ESP1010FNR2		10		
ESP1012FNR2		12		
ESP1014FNR2		14		
ESP1015FNR2		15		
ESP1206FNR2	12	6	14	70
ESP1208FNR2		8		
ESP1210FNR2		10		
ESP1212FNR2		12		
ESP1214FNR2		14		
ESP1215FNR2		15		
ESP1406FNR2	14	6	14	70
ESP1408FNR2		8		
ESP1410FNR2		10		
ESP1412FNR2		12		
ESP1414FNR2		14		
ESP1415FNR2		15		
ESP1606FNR2	16	6	16	70
ESP1608FNR2		8		
ESP1610FNR2		10		
ESP1612FNR2		12		
ESP1614FNR2		14		
ESP1615FNR2		15		
ESP1806FNR2	18	6	16	70
ESP1808FNR2		8		
ESP1810FNR2		10		
ESP1812FNR2		12		
ESP1814FNR2		14		
ESP1815FNR2		15		

Proximal Release				
Distal & Proximal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ESP2006FNR2	20	6	20	70
ESP2008FNR2		8		
ESP2010FNR2		10		
ESP2012FNR2		12		
ESP2014FNR2		14		
ESP2015FNR2		15		
ESP2206FNR2	22	6	20	70
ESP2208FNR2		8		
ESP2210FNR2		10		
ESP2212FNR2		12		
ESP2214FNR2		14		
ESP2215FNR2		15		
ESP2406FNR2	24	6	22	70
ESP2408FNR2		8		
ESP2410FNR2		10		
ESP2412FNR2		12		
ESP2414FNR2		14		
ESP2415FNR2		15		
ESP2806FNR2	28	6	22	70
ESP2808FNR2		8		
ESP2810FNR2		10		
ESP2812FNR2		12		
ESP2814FNR2		14		
ESP2815FNR2		15		