


# UVENTA™ Urethral Stent

for maintaining urethral patency in benign and malignant urethral obstruction



## Features

- Fully-covered with Silicone
- Self-expandable
- Anti-migration barbs
- Radiopaque markers
- Retrieval Suture

# Ordering Information

TTS (Through The Scope)							
Straight							
Code	Stent		Delivery				
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)			
US1003F	10	3	10	50			
US1004F		4					
US1005F		5					
US1006F		6					
US1007F		7					
US1008F		8					
US1009F		9					
US1010F		10					
US1203F		12			3	10	50
US1204F					4		
US1205F	5						
US1206F	6						
US1207F	7						
US1208F	8						
US1209F	9						
US1210F	10						

TTS (Through The Scope)							
Straight							
Code	Stent		Delivery				
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)			
US1403F	14	3	10	50			
US1404F		4					
US1405F		5					
US1406F		6					
US1407F		7					
US1408F		8					
US1409F		9					
US1410F		10					
US1603F		16			3	10	50
US1604F					4		
US1605F	5						
US1606F	6						
US1607F	7						
US1608F	8						
US1609F	9						
US1610F	10						