

Niti-S Esophageal Covered Stent [Both Bare-Type]



Niti-S Esophageal Covered Stent [Both Bare-Type]

for malignant esophageal strictures

Ordering Information

Distal Release				
Both ends bare, Distal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ES1006B	10	6	12	70
ES1008B		8		
ES1010B		10		
ES1012B		12		
ES1014B		14		
ES1015B		15		
ES1206B	12	6	14	70
ES1208B		8		
ES1210B		10		
ES1212B		12		
ES1214B		14		
ES1215B		15		
ES1406B	14	6	14	70
ES1408B		8		
ES1410B		10		
ES1412B		12		
ES1414B		14		
ES1415B		15		
ES1606B	16	6	16	70
ES1608B		8		
ES1610B		10		
ES1612B		12		
ES1614B		14		
ES1615B		15		
ES1806B	18	6	16	70
ES1808B		8		
ES1810B		10		
ES1812B		12		
ES1814B		14		
ES1815B		15		

Distal Release				
Both ends bare, Distal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ES2006B	20	6	20	70
ES2008B		8		
ES2010B		10		
ES2012B		12		
ES2014B		14		
ES2015B		15		
ES2206B	22	6	20	70
ES2208B		8		
ES2210B		10		
ES2212B		12		
ES2214B		14		
ES2215B		15		
ES2406B	24	6	22	70
ES2408B		8		
ES2410B		10		
ES2412B		12		
ES2414B		14		
ES2415B		15		
ES2806B	28	6	22	70
ES2808B		8		
ES2810B		10		
ES2812B		12		
ES2814B		14		
ES2815B		15		

Ordering Information

Proximal Release				
Both ends bare, Distal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ESP1006B	10	6	14	70
ESP1008B		8		
ESP1010B		10		
ESP1012B		12		
ESP1014B		14		
ESP1015B		15		
ESP1206B	12	6	14	70
ESP1208B		8		
ESP1210B		10		
ESP1212B		12		
ESP1214B		14		
ESP1215B		15		
ESP1406B	14	6	14	70
ESP1408B		8		
ESP1410B		10		
ESP1412B		12		
ESP1414B		14		
ESP1415B		15		
ESP1606B	16	6	16	70
ESP1608B		8		
ESP1610B		10		
ESP1612B		12		
ESP1614B		14		
ESP1615B		15		
ESP1806B	18	6	16	70
ESP1808B		8		
ESP1810B		10		
ESP1812B		12		
ESP1814B		14		
ESP1815B		15		

Proximal Release				
Both ends bare, Distal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
ESP2006B	20	6	20	70
ESP2008B		8		
ESP2010B		10		
ESP2012B		12		
ESP2014B		14		
ESP2015B		15		
ESP2206B	22	6	20	70
ESP2208B		8		
ESP2210B		10		
ESP2212B		12		
ESP2214B		14		
ESP2215B		15		
ESP2406B	24	6	22	70
ESP2408B		8		
ESP2410B		10		
ESP2412B		12		
ESP2414B		14		
ESP2415B		15		
ESP2806B	28	6	22	70
ESP2808B		8		
ESP2810B		10		
ESP2812B		12		
ESP2814B		14		
ESP2815B		15		

Ordering Information

TTS (Through The Scope)				
Both ends bare, Distal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
EST1806B	18	6	10.5	160
EST1808B		8		
EST1810B		10		
EST1812B		12		
EST1814B		14		
EST1815B		15		
EST2006B	20	6	10.5	160
EST2008B		8		
EST2010B		10		
EST2012B		12		
EST2014B		14		
EST2015B		15		
EST2206B	22	6	10.5	160
EST2208B		8		
EST2210B		10		
EST2212B		12		
EST2214B		14		
EST2215B		15		

TTS (Through The Scope)				
Both ends bare, Distal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
EST1806B-18	18	6	10.5	180
EST1808B-18		8		
EST1810B-18		10		
EST1812B-18		12		
EST1814B-18		14		
EST1815B-18		15		
EST2006B-18	20	6	10.5	180
EST2008B-18		8		
EST2010B-18		10		
EST2012B-18		12		
EST2014B-18		14		
EST2015B-18		15		
EST2206B-18	22	6	10.5	180
EST2208B-18		8		
EST2210B-18		10		
EST2212B-18		12		
EST2214B-18		14		
EST2215B-18		15		

TTS (Through The Scope)				
Both ends bare, Distal Removal String				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
EST1806B-22	18	6	10.5	220
EST1808B-22		8		
EST1810B-22		10		
EST1812B-22		12		
EST1814B-22		14		
EST1815B-22		15		
EST2006B-22	20	6	10.5	220
EST2008B-22		8		
EST2010B-22		10		
EST2012B-22		12		
EST2014B-22		14		
EST2015B-22		15		
EST2206B-22	22	6	10.5	220
EST2208B-22		8		
EST2210B-22		10		
EST2212B-22		12		
EST2214B-22		14		
EST2215B-22		15		