


Niti-S Pyloric/Duodenal Uncovered Stent [D-Type]

for intrinsic and/or extrinsic malignant pyloric/duodenal stricture



Ordering Information

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PD1806	18	6	12	135
PD1808		8		
PD1810		10		
PD1812		12		
PD1814		14		
PD1815		15		
PD2006	20	6	12	135
PD2008		8		
PD2010		10		
PD2012		12		
PD2014		14		
PD2015		15		
PD2206	22	6	16	135
PD2208		8		
PD2210		10		
PD2212		12		
PD2214		14		
PD2215		15		
PD2406	24	6	16	135
PD2408		8		
PD2410		10		
PD2412		12		
PD2414		14		
PD2415		15		
PD2606	26	6	18	135
PD2608		8		
PD2610		10		
PD2612		12		
PD2614		14		
PD2615		15		
PD2806	28	6	18	135
PD2808		8		
PD2810		10		
PD2812		12		
PD2814		14		
PD2815		15		

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PDT1806	18	6	10	180
PDT1808		8		
PDT1810		10		
PDT1812		12		
PDT1814		14		
PDT1815		15		
PDT2006	20	6	10	180
PDT2008		8		
PDT2010		10		
PDT2012		12		
PDT2014		14		
PDT2015		15		
PDT2206	22	6	10	180
PDT2208		8		
PDT2210		10		
PDT2212		12		
PDT2214		14		
PDT2215		15		
PDT2406	24	6	10	180
PDT2408		8		
PDT2410		10		
PDT2412		12		
PDT2414		14		
PDT2415		15		

Ordering Information

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PDT1806-12	18	6	10	120
PDT1808-12		8		
PDT1810-12		10		
PDT1812-12		12		
PDT1814-12		14		
PDT1815-12		15		
PDT2006-12	20	6	10	120
PDT2008-12		8		
PDT2010-12		10		
PDT2012-12		12		
PDT2014-12		14		
PDT2015-12		15		
PDT2206-12	22	6	10	120
PDT2208-12		8		
PDT2210-12		10		
PDT2212-12		12		
PDT2214-12		14		
PDT2215-12		15		
PDT2406-12	24	6	10	120
PDT2408-12		8		
PDT2410-12		10		
PDT2412-12		12		
PDT2414-12		14		
PDT2415-12		15		

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PDT1806-22	18	6	10	220
PDT1808-22		8		
PDT1810-22		10		
PDT1812-22		12		
PDT1814-22		14		
PDT1815-22		15		
PDT2006-22	20	6	10	220
PDT2008-22		8		
PDT2010-22		10		
PDT2012-22		12		
PDT2014-22		14		
PDT2015-22		15		
PDT2206-22	22	6	10	220
PDT2208-22		8		
PDT2210-22		10		
PDT2212-22		12		
PDT2214-22		14		
PDT2215-22		15		
PDT2406-22	24	6	10	220
PDT2408-22		8		
PDT2410-22		10		
PDT2412-22		12		
PDT2414-22		14		
PDT2415-22		15		

Ordering Information

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PDT1806-F9	18	6	9	180
PDT1808-F9		8		
PDT1810-F9		10		
PDT1812-F9		12		
PDT1814-F9		14		
PDT1815-F9		15		
PDT2006-F9	20	6	9	180
PDT2008-F9		8		
PDT2010-F9		10		
PDT2012-F9		12		
PDT2014-F9		14		
PDT2015-F9		15		
PDT2206-F9	22	6	9	180
PDT2208-F9		8		
PDT2210-F9		10		
PDT2212-F9		12		
PDT2214-F9		14		
PDT2215-F9		15		
PDT2406-F9	24	6	9	180
PDT2408-F9		8		
PDT2410-F9		10		
PDT2412-F9		12		
PDT2414-F9		14		
PDT2415-F9		15		

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PDT1806-22F9	18	6	9	220
PDT1808-22F9		8		
PDT1810-22F9		10		
PDT1812-22F9		12		
PDT1814-22F9		14		
PDT1815-22F9		15		
PDT2006-22F9	20	6	9	220
PDT2008-22F9		8		
PDT2010-22F9		10		
PDT2012-22F9		12		
PDT2014-22F9		14		
PDT2015-22F9		15		
PDT2206-22F9	22	6	9	220
PDT2208-22F9		8		
PDT2210-22F9		10		
PDT2212-22F9		12		
PDT2214-22F9		14		
PDT2215-22F9		15		
PDT2406-22F9	24	6	9	220
PDT2408-22F9		8		
PDT2410-22F9		10		
PDT2412-22F9		12		
PDT2414-22F9		14		
PDT2415-22F9		15		