

# Niti-S Enteral Colonic Covered Stent [Both Bare-Type]



# Niti-S Enteral Colonic Covered Stent [Both Bare-Type]

for intrinsic and/or extrinsic malignant colonic stricture

## Ordering Information

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CST1806B	18	6	10.5	220
CST1808B		8		
CST1810B		10		
CST1812B		12		
CST1814B		14		
CST1815B		15		
CST2006B	20	6	10.5	220
CST2008B		8		
CST2010B		10		
CST2012B		12		
CST2014B		14		
CST2015B		15		
CST2206B	22	6	10.5	220
CST2208B		8		
CST2210B		10		
CST2212B		12		
CST2214B		14		
CST2215B		15		

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CS2006B	20	6	20	70
CS2008B		8		
CS2010B		10		
CS2012B		12		
CS2014B		14		
CS2015B		15		
CS2016B	22	16	20	70
CS2206B		6		
CS2208B		8		
CS2210B		10		
CS2212B		12		
CS2214B		14		
CS2215B	24	15	22	70
CS2216B		16		
CS2406B		6		
CS2408B		8		
CS2410B		10		
CS2412B		12		
CS2414B	26	14	22	70
CS2415B		15		
CS2416B		16		
CS2606B		6		
CS2608B		8		
CS2610B		10		
CS2612B	28	12	22	70
CS2614B		14		
CS2615B		15		
CS2616B		16		
CS2806B		6		
CS2808B		8		
CS2810B	28	10	22	70
CS2812B		12		
CS2814B		14		
CS2815B		15		
CS2816B		16		

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
CS1806B	18	6	16	70
CS1808B		8		
CS1810B		10		
CS1812B		12		
CS1814B		14		
CS1815B		15		
CS1816B		16		