


# Niti-S Nagi™ Stent

For drainage of a pancreatic pseudocyst through a transgastric or transduodenal approach



## Features

- **Wide and smooth flare edges**
  - Prevent the risk of migration
- **Available in various diameters (Up to 16mm)**
- **Retrieval String** for repositioning or easy removal
- **Radiopaque marker:** 3 (three) at both ends & 2 (two) in the middle

# Ordering Information

Endoscopic Approach				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BS1001FW	10	1	9	180
BS1002FW		2		
BS1003FW		3		
BS1201FW	12	1	9	180
BS1202FW		2		
BS1203FW		3		
BS1402FW	14	2	10	180
BS1403FW		3		
BS1602FW	16	2	10	180
BS1603FW		3		