


Niti-S Biliary Uncovered Stent [LCD-Type]

for malignant biliary strictures



Features

- **Unfixed cell with a weaving construction**
 - Low axial force to maintain luminal patency in tortuous anatomy
- **Radiopaque marker: 3 (three) at both ends & 2 (two) in the middle**

Ordering Information

Percutaneous				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TLD0604	6	4	8	50
TLD0605		5		
TLD0606		6		
TLD0607		7		
TLD0608		8		
TLD0609		9		
TLD0610		10		
TLD0612		12		
TLD0804	8	4	8	50
TLD0805		5		
TLD0806		6		
TLD0807		7		
TLD0808		8		
TLD0809		9		
TLD0810		10		
TLD0812		12		
TLD1004	10	4	8	50
TLD1005		5		
TLD1006		6		
TLD1007		7		
TLD1008		8		
TLD1009		9		
TLD1010		10		
TLD1012		12		

Percutaneous				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TLD0604-6	6	4	6	50
TLD0605-6		5		
TLD0606-6		6		
TLD0607-6		7		
TLD0608-6		8		
TLD0609-6		9		
TLD0610-6		10		
TLD0612-6		12		
TLD0804-6	8	4	6	50
TLD0805-6		5		
TLD0806-6		6		
TLD0807-6		7		
TLD0808-6		8		
TLD0809-6		9		
TLD0810-6		10		
TLD0812-6		12		

Percutaneous				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TLD0604-7	6	4	7	50
TLD0605-7		5		
TLD0606-7		6		
TLD0607-7		7		
TLD0608-7		8		
TLD0609-7		9		
TLD0610-7		10		
TLD0612-7		12		
TLD0804-7	8	4	7	50
TLD0805-7		5		
TLD0806-7		6		
TLD0807-7		7		
TLD0808-7		8		
TLD0809-7		9		
TLD0810-7		10		
TLD0812-7		12		
TLD1004-7	10	4	7	50
TLD1005-7		5		
TLD1006-7		6		
TLD1007-7		7		
TLD1008-7		8		
TLD1009-7		9		
TLD1010-7		10		
TLD1012-7		12		

Percutaneous				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TLD1004-6	10	4	6	50
TLD1005-6		5		
TLD1006-6		6		
TLD1007-6		7		
TLD1008-6		8		
TLD1009-6		9		
TLD1010-6		10		
TLD1012-6		12		

Ordering Information

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BLD0604	6	4	8	180
BLD0605		5		
BLD0606		6		
BLD0607		7		
BLD0608		8		
BLD0609		9		
BLD0610		10		
BLD0612		12		
BLD0804	8	4	8	180
BLD0805		5		
BLD0806		6		
BLD0807		7		
BLD0808		8		
BLD0809		9		
BLD0810		10		
BLD0812		12		
BLD1004	10	4	8	180
BLD1005		5		
BLD1006		6		
BLD1007		7		
BLD1008		8		
BLD1009		9		
BLD1010		10		
BLD1012		12		

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BLD0604-7	6	4	7	180
BLD0605-7		5		
BLD0606-7		6		
BLD0607-7		7		
BLD0608-7		8		
BLD0609-7		9		
BLD0610-7		10		
BLD0612-7		12		
BLD0804-7	8	4	7	180
BLD0805-7		5		
BLD0806-7		6		
BLD0807-7		7		
BLD0808-7		8		
BLD0809-7		9		
BLD0810-7		10		
BLD0812-7		12		
BLD1004-7	10	4	7	180
BLD1005-7		5		
BLD1006-7		6		
BLD1007-7		7		
BLD1008-7		8		
BLD1009-7		9		
BLD1010-7		10		
BLD1012-7		12		

Ordering Information

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BLD0604-6	6	4	6	180
BLD0605-6		5		
BLD0606-6		6		
BLD0607-6		7		
BLD0608-6		8		
BLD0609-6		9		
BLD0610-6		10		
BLD0612-6		12		
BLD0804-6	8	4	6	180
BLD0805-6		5		
BLD0806-6		6		
BLD0807-6		7		
BLD0808-6		8		
BLD0809-6		9		
BLD0810-6		10		
BLD0812-6		12		
BLD1004-6	10	4	6	180
BLD1005-6		5		
BLD1006-6		6		
BLD1007-6		7		
BLD1008-6		8		
BLD1009-6		9		
BLD1010-6		10		
BLD0604-6		12		

Endoscopic, Short-wire				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BLDM0604	6	4	8.5	180
BLDM0605		5		
BLDM0606		6		
BLDM0607		7		
BLDM0608		8		
BLDM0609		9		
BLDM0610		10		
BLDM0804	8	4	8.5	180
BLDM0805		5		
BLDM0806		6		
BLDM0807		7		
BLDM0808		8		
BLDM0809		9		
BLDM0810		10		
BLDM1004	10	4	8.5	180
BLDM1005		5		
BLDM1006		6		
BLDM1007		7		
BLDM1008		8		
BLDM1009		9		
BLDM1010		10		