


Niti-S Pyloric/Duodenal Covered Stent



Features

- **Fixed cell with braided construction**
- **Silicone covering and soft round ends**
 - Reduce tissue ingrowth
- **Visible green suture for easy removal**
- **Radiopaque marker: 4 (four) at both ends & 2 (two) in the middle**

Niti-S Pyloric/Duodenal Covered Stent [Full Covered-Type]

for intrinsic and/or extrinsic benign and/or malignant pyloric/duodenal stricture

Ordering Information

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PS1806F	18	6	16	135
PS1808F		8		
PS1810F		10		
PS1812F		12		
PS1814F		14		
PS1815F		15		
PS1816F		16		
PS2006F	20	6	16	135
PS2008F		8		
PS2010F		10		
PS2012F		12		
PS2014F		14		
PS2015F		15		
PS2016F		16		
PS2206F	22	6	20	135
PS2208F		8		
PS2210F		10		
PS2212F		12		
PS2214F		14		
PS2215F		15		
PS2216F		16		

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PS2406F	24	6	22	135
PS2408F		8		
PS2410F		10		
PS2412F		12		
PS2414F		14		
PS2415F		15		
PS2416F		16		
PS2606F	26	6	22	135
PS2608F		8		
PS2610F		10		
PS2612F		12		
PS2614F		14		
PS2615F		15		
PS2616F		16		
PS2806F	28	6	22	135
PS2808F		8		
PS2810F		10		
PS2812F		12		
PS2814F		14		
PS2815F		15		
PS2816F		16		

Ordering Information

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PST1806F	18	6	10.5	180
PST1808F		8		
PST1810F		10		
PST1812F		12		
PST1814F		14		
PST1815F		15		
PST2006F	20	6	10.5	180
PST2008F		8		
PST2009F		9		
PST2010F		10		
PST2012F		12		
PST2014F		14		
PST2015F		15		
PST2206F	22	6	10.5	180
PST2208F		8		
PST2209F		9		
PST2210F		10		
PST2212F		12		
PST2214F		14		
PST2215F		15		

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PST1806F-12	18	6	10.5	120
PST1808F-12		8		
PST1810F-12		10		
PST1812F-12		12		
PST1814F-12		14		
PST1815F-12		15		
PST2006F-12	20	6	10.5	120
PST2008F-12		8		
PST2010F-12		10		
PST2012F-12		12		
PST2014F-12		14		
PST2015F-12		15		
PST2206F-12	22	6	10.5	120
PST2208F-12		8		
PST2210F-12		10		
PST2212F-12		12		
PST2214F-12		14		
PST2215F-12		15		

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PST1806F-22	18	6	10.5	220
PST1808F-22		8		
PST1810F-22		10		
PST1812F-22		12		
PST1814F-22		14		
PST1815F-22		15		
PST2006F-22	20	6	10.5	220
PST2008F-22		8		
PST2010F-22		10		
PST2012F-22		12		
PST2014F-22		14		
PST2015F-22		15		
PST2206F-22	22	6	10.5	220
PST2208F-22		8		
PST2210F-22		10		
PST2212F-22		12		
PST2214F-22		14		
PST2215F-22		15		

Niti-S Pyloric/Duodenal Covered Stent [Both Bare-Type]

for intrinsic and/or extrinsic malignant pyloric/duodenal stricture

Ordering Information

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PS1806B	18	6	16	135
PS1808B		8		
PS1810B		10		
PS1812B		12		
PS1814B		14		
PS1815B		15		
PS1816B		16		
PS2006B	20	6	16	135
PS2008B		8		
PS2010B		10		
PS2012B		12		
PS2014B		14		
PS2015B		15		
PS2016B		16		
PS2206B	22	6	20	135
PS2208B		8		
PS2210B		10		
PS2212B		12		
PS2214B		14		
PS2215B		15		
PS2216B		16		

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PS2406B	24	6	22	135
PS2408B		8		
PS2410B		10		
PS2412B		12		
PS2414B		14		
PS2415B		15		
PS2416B		16		
PS2606B	26	6	22	135
PS2608B		8		
PS2610B		10		
PS2612B		12		
PS2614B		14		
PS2615B		15		
PS2616B		16		
PS2806B	28	6	22	135
PS2808B		8		
PS2810B		10		
PS2812B		12		
PS2814B		14		
PS2815B		15		
PS2816B		16		

Ordering Information

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PST1806B	18	6	10.5	180
PST1808B		8		
PST1810B		10		
PST1812B		12		
PST1814B		14		
PST1815B		15		
PST2006B	20	6	10.5	180
PST2008B		8		
PST2009B		9		
PST2010B		10		
PST2012B		12		
PST2014B		14		
PST2015B		15		
PST2206B	22	6	10.5	180
PST2208B		8		
PST2209B		9		
PST2210B		10		
PST2212B		12		
PST2214B		14		
PST2215B		15		

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PST1806B-12	18	6	10.5	120
PST1808B-12		8		
PST1810B-12		10		
PST1812B-12		12		
PST1814B-12		14		
PST1815B-12		15		
PST2006B-12	20	6	10.5	120
PST2008B-12		8		
PST2010B-12		10		
PST2012B-12		12		
PST2014B-12		14		
PST2015B-12		15		
PST2206B-12	22	6	10.5	120
PST2208B-12		8		
PST2210B-12		10		
PST2212B-12		12		
PST2214B-12		14		
PST2215B-12		15		

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PST1806B-22	18	6	10.5	220
PST1808B-22		8		
PST1810B-22		10		
PST1812B-22		12		
PST1814B-22		14		
PST1815B-22		15		
PST2006B-22	20	6	10.5	220
PST2008B-22		8		
PST2010B-22		10		
PST2012B-22		12		
PST2014B-22		14		
PST2015B-22		15		
PST2206B-22	22	6	10.5	220
PST2208B-22		8		
PST2210B-22		10		
PST2212B-22		12		
PST2214B-22		14		
PST2215B-22		15		

Niti-S Pyloric/Duodenal Covered Stent [End Bare-Type]

for intrinsic and/or extrinsic malignant pyloric/duodenal stricture

Ordering Information

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PS1806H	18	6	16	135
PS1808H		8		
PS1810H		10		
PS1812H		12		
PS1814H		14		
PS1815H		15		
PS1816H		16		
PS2006H	20	6	16	135
PS2008H		8		
PS2010H		10		
PS2012H		12		
PS2014H		14		
PS2015H		15		
PS2016H		16		
PS2206H	22	6	20	135
PS2208H		8		
PS2210H		10		
PS2212H		12		
PS2214H		14		
PS2215H		15		
PS2216H		16		

OTW (Over The Wire)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PS2406H	24	6	22	135
PS2408H		8		
PS2410H		10		
PS2412H		12		
PS2414H		14		
PS2415H		15		
PS2416H		16		
PS2606H	26	6	22	135
PS2608H		8		
PS2610H		10		
PS2612H		12		
PS2614H		14		
PS2615H		15		
PS2616H		16		
PS2806H	28	6	22	135
PS2808H		8		
PS2810H		10		
PS2812H		12		
PS2814H		14		
PS2815H		15		
PS2816H		16		

Ordering Information

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PST1806H	18	6	10.5	180
PST1808H		8		
PST1810H		10		
PST1812H		12		
PST1814H		14		
PST1815H		15		
PST2006H	20	6	10.5	180
PST2008H		8		
PST2010H		10		
PST2012H		12		
PST2014H		14		
PST2015H		15		
PST2206H	22	6	10.5	180
PST2208H		8		
PST2210H		10		
PST2212H		12		
PST2214H		14		
PST2215H		15		

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PST1806H-12	18	6	10.5	120
PST1808H-12		8		
PST1810H-12		10		
PST1812H-12		12		
PST1814H-12		14		
PST1815H-12		15		
PST2006H-12	20	6	10.5	120
PST2008H-12		8		
PST2010H-12		10		
PST2012H-12		12		
PST2014H-12		14		
PST2015H-12		15		
PST2206H-12	22	6	10.5	120
PST2208H-12		8		
PST2210H-12		10		
PST2212H-12		12		
PST2214H-12		14		
PST2215H-12		15		

TTS (Through The Scope)				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
PST1806H-22	18	6	10.5	220
PST1808H-22		8		
PST1810H-22		10		
PST1812H-22		12		
PST1814H-22		14		
PST1815H-22		15		
PST2006H-22	20	6	10.5	220
PST2008H-22		8		
PST2010H-22		10		
PST2012H-22		12		
PST2014H-22		14		
PST2015H-22		15		
PST2206H-22	22	6	10.5	220
PST2208H-22		8		
PST2210H-22		10		
PST2212H-22		12		
PST2214H-22		14		
PST2215H-22		15		