

Optimos™

Cystotome

is intended to dilate transmural punctures to facilitate endoscopic access to the gastrointestinal lumen or tract, or cystic cavities



Product No.	Sheath O.D.	Tip O.D.	Usable Length	Recommended Guide wire
OCT1906	Ø2.0mm / 6Fr	Ø1.4mm	1,850mm	0.025~0.035"
OCT1907	Ø2.3mm / 7Fr	Ø1.4mm		
OCT1908	Ø2.8mm / 8.5Fr	Ø1.7mm		

HEAD OFFICE
14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, Korea 10022

SALES OFFICE
Bldg.C 4F, 158 Haneulmaeul-ro, Ilsandong-gu, Goyang-si, Gyeonggi-do, Korea 10355
E: contact@stent.net T: +82 31 904 6153 F: +82 31 904 6157

www.taewoongmedical.com

