


Niti-S Biliary Covered Stent [Giobor]

for malignant biliary strictures



Features

- **Partially covered design**
 - (1) The long covered portion (70%)
 - (2) The bare portion (30%)
- **Anti-migration flared end**
- **Radiopaque markers:** Three (3) at both ends & Two (2) at the end of the covered part

Ordering Information

Endoscopic				
Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
BS0806FP	8	6	8.5	180
BS0808FP		8		
BS0810FP		10		
BS1006FP	10	6	8.5	180
BS1008FP		8		
BS1010FP		10		